

# Protocol for probing DS

7th May 2024

# Stimuli

- Sweeping bar (ON/OFF) (V/H) (1pix / 4pix)
- Sweeping edge (ON/OFF) (V/H)
- Sweeping diagonal bar (ON/OFF) (D/-D) (4pix)
  
- Full field flash (ON/OFF)
- Gratings (ON/OFF) (V/H/D/-D)

*N.B. both directions*

# Stimuli

- Sweeping bar (ON/OFF) (V/H) (1pix / 4pix)
  - Sweeping edge (ON/OFF) (V/H)
  - Sweeping diagonal bar (ON/OFF) (D/-D) (4pix)
  - Full field flash (ON/OFF)
  - Gratings (ON/OFF) (V/H/D/-D)
- 20 deg/s, 100 deg/s
  - 20 deg/s, 100 deg/s
  - 20 deg/s, 100 deg/s
  - 0.5Hz, 1Hz, 2Hz, 4Hz, 8Hz, 16Hz, 32HZ, 64Hz, 75Hz
  - 0.5Hz, 1Hz, 2Hz, 4Hz, 8Hz, 16Hz, 32HZ, 64Hz, 75Hz




*N.B. both directions*

# Stimuli

- Sweeping bar (ON/OFF) (V/H) (1pix / 4pix)
- Sweeping edge (ON/OFF) (V/H)
- Sweeping diagonal bar (ON/OFF) (D/-D) (4pix)

- 20 deg/s, 100 deg/s
- 20 deg/s, 100 deg/s
- 20 deg/s, 100 deg/s











*Isaacson et al. 2022*

- 0.5Hz, 1Hz, 2Hz, 4Hz, 8Hz, 16Hz, 32Hz, 64Hz, 75Hz
- 0.5Hz, 1Hz, 2Hz, 4Hz, 8Hz, 16Hz, 32Hz, 64Hz, 75Hz



- Sweeping bar (ON/OFF) (V/H) (1pix / 4pix)




20: **3s**, 100:**0.8s**

20: **7s**, 100:**1.6s**




- Sweeping edge (ON/OFF) (V/H)



20: **2.8s**, 100:**0.6s**







20: **7s**, 100:**1.6s**









- Full field flash (ON/OFF)



**3s** per speed - 9 speeds = **27s** each







- Gratings (ON/OFF) (V/H)



**3.25s** per speed - 9 speeds = **30s** each






- Sweeping diagonal bar (ON/OFF) (D/-D) (4pix)



20: **5.1s**, 100: **1.2s** - need to check



- Gratings (ON/OFF) (D/-D)



**3.25s** per speed - 9 speeds = **30s** each



- 0.25s before beginning of stimuli (gratings, bars)
  - 36 minutes for 3 rep
- ~12 mins per rep
- Using 0, 7, 14 for light intensity - too bright?
  - 1 pix diagonal bar?
  - Sweeping stimuli only cover 135 deg.

